



**MEDICAL PRODUCT CATALOGUE**

# OFFICIAL PHILIPPINE PARTNER TO



**AND**  
Tokyo, Japan

Sphygmomanometer &  
Digital BP Devices



**Ambu**  
Denmark

Video Laryngoscope



**CHINESPORT**  
ITALY

Italy  
Rehabilitation Equipment /  
Therapy Tables /  
Patient Transport



**NAVEH**  
PHARMA

Israel  
Earwax Removal Spray



**OPHARDT**  
HYGIENE

Germany  
Hand Hygiene  
Dispensers



**seca**  
Precision for health

Hamburg, Germany  
Precise Medical  
Measurements &  
Weighing Scales



**G/NEVR/**  
Neo-Tech for Life

Italy  
Incubators and  
Neonate Device



**Guldmann**  
DENMARK

Patient Lifters



**Hadeco**  
Kawasaki, Japan

Fetal, Vascular & Surgical  
Dopplers



**SECULIFE**

Germany  
Medical Device Testers  
/ Simulators



**Schaerer**  
MEDICAL

Switzerland  
Mobile Operating Tables



**SINAPI**  
biomedical

South Africa  
Chest Drainage System



**KaWe**  
GERMANY

Laryngoscope  
Otoscope  
Ophthalmoscope  
Diagnostic Sets



**KOIKE MEDICAL**

Japan  
User-Friendly  
Oxygen Regulator



**MAICO**

Germany  
Newborn Hearing  
Screening Devices  
OAE/AABR



**SURGIRIS**  
SMART VISION

France  
Surgical Lights



**TECNIMED**

Italy  
Touchless Thermometer



**TRISMED**  
Cardiac Diagnosis Systems

South Korea  
Cardiac Diagnosis Systems

# GREETINGS!

We bring you the latest publication of the Wellness PRO Product Book with much excitement and anticipation. This book has been put together in the most interesting of times to date.

We formulated this book with COVID-19 (iterations, surges, variants, and all) raging across the globe, including the Philippines.

The pandemic has been especially poignant for the Philippines. We garnered the distinction of having the most extended lockdown in the world. Covid-19 forced a reckoning and a measuring of our Healthcare system, and we have, unfortunately, been found wanting.

We could not even take care of our healers and healthcare professionals when Covid stricken initially. As we locked down in March 2020, we found only around 2,000 ICU beds for a population of about 109 million people. We had to review hand washing that we already knew since kindergarten. We were chased into our homes – faces covered, family meals avoided, coughs, colds and chills challenged respirations, and occasional hypoxia nursed. Meanwhile, we got educated on the Greek alphabet, spike proteins, mRNAs, comorbidities, RT-PCR, immunity, and ventilation as a nation.

One clear message of this pandemic to the Philippine Healthcare is “Better Shape Up!” be it in terms of availability of medical personnel, infrastructure, information exchange, funding, patient education, political will, of medical statistics, among others.

It is in response to the above realization that we put forth this Wellness PRO Product book.

We are of the belief that there are indeed medical practices and medical technology solutions, proven as better in terms of quality patient outcome and already utilized in a mainstream way, around the world. These must be made available to our Filipino patients, too, sooner than later. Truth be told, to a large extent, access to these is conditioned only in part by the idea of funding, but more by the (limited) extent we view Philippine healthcare strategically vis a vis our vision for our nation.

The pandemic should come to pass, sooner or later. Our drive to improve Philippine Healthcare is a work in progress that we all have to undertake today and for generations to come.

We hope that you find this Wellness PRO Product Book enjoyable and valuable.

Take Care! Stay Safe! Stay Well!



**JOSEPH CHAN JUAN**  
PRESIDENT



## KingVISION® aBlade Video Laryngoscope

### ✦ PORTABLE

- King Vision is light weight, self-contained, battery operated and water resistant.
- King Vision can be used continuously for more than 90 minutes making it highly portable.
- Blades are individually packaged. Each blade is disposable, eliminating concerns about cross contamination.
- King Vision has a High Intensity OLED display that creates clear image viewing in 160° panoramic field which is also beneficial in teaching situations.

### ✦ AFFORDABLE

- Disposable aBlades allow economical use of the King Vision video laryngoscope for all intubations. With the addition of the aBlade configuration you have an excellent solution to the problem of limited availability of advanced imaging devices in the OR.
- With the King Vision you can improve patient care by having high performance visualisation capabilities for every intubation at an affordable price.

### ✦ DURABLE

- Designed to be your primary tool for all intubations.
- Video output capability is compatible with external monitors and recording devices.
- The robust, full-colour, non-glare display can resist repeated cleaning and normal use wear and tear. The video adapter is made from durable materials that can withstand repeated use.



## Key Benefits

Experience the possibilities



### High flexibility

King Vision offers the choice of two blade types for user flexibility.



### Operational in seconds

King Vision is light weight, self-contained and battery operated, which makes it operational in seconds.



### Minimal impact

King Vision is ergonomically designed to provide minimal lifting of soft tissue and impact on teeth.



### Visionary performance

King Vision holds a full colour TFT LCD display with advanced imaging technology that provides visionary performance in the palm of your hand.



### No cross-contamination

Single-use blades eliminate risk of cross contamination.




### Training

With a clear image viewing in 160° panoramic field King Vision is suitable for training of intubation skills.

## KingVISION® aBlade Video Laryngoscope

### Technical Specification

REUSABLE DISPLAY	
Batteries	3 AAA (Alkaline recommended)
Battery Life	> 90 min. (Battery status indicator flashes red when batteries need to be changed.)
Computerized Power Management System	Auto shut off, Auto white balancing
DISPOSABLE BLADES	
Anti-fog	Coating on distal lens
AP Height of Blade	13 mm (standard); 18 mm (channeled)
Blade Length	17 cm
Blade Width	26 mm (standard); 29 mm (channeled); 16 mm at distal tip
Disposable Material	Polycarbonate / TPE
ET Size	6.0 mm - 8.0 mm (Channeled Blade)
REUSABLE VIDEO ADAPTER	
Video Aspect Ratio	4:3
Video Port	RCA connection to monitor with cable
Video Refresh Rate	30 frames per second
Video Resolution	320 x 240 (QVGA)
Video Screen	OLED Display
Video Screen Size	6.1 cm / 2.4" diagonal
Material	Polycarbonate / ABS anti-reflective coating on display window
Camera Chip	CMOS
Camera Resolution	640 x 480 VGA
Light Source	White LED




The Universal Table for All Surgical Discipline

# AXIS 500

## MOBILE OPERATING TABLE

The Universal Table



The Schaerer AXIS 500 is the most successful product line of Schaerer.

The universal OT-table is an universal OT-table for all surgical disciplines with one electro-hydraulically driven joint. Removable two-part seat plate bars to connect high specialized accessories.

Due to the unique cantilevered column design AXIS offers the largest X-ray window of 1450mm / 57" with 220kg / 485lbs patient weight capacity and unobstructed C-arm access. There is no sliding needed. AXIS is the ideal universal OT-imaging table and set the standard is this range.

## Key Features

- Ergonomic remote control panel
- Emergency built-in control panel
- Visco elastic pods (reduces pressure up to 50% compared to standard)
- Electro-hydraulically driven table top
- Removable, two-part seat plate bars for fitting orthopedic accessories
- Manually adjustable head and leg plates supported by gas pressure springs
- Universal for all accessories

## Maximum quality with Swiss precision

Manufactured in Switzerland with design and engineering on highest standards, in accordance with certified quality standards ISO 13485:2016.

## Unique Benefits of AXIS Table

### MAXIMUM STABILITY

The use of high-quality Swiss materials and components ensures strength and a long service life. 220 kg / 485 lbs all position dynamic weight capacity.

### MAXIMUM SAFETY

Only Schaerer operating tables have two independent motor pump systems that can work powered by batteries for several operating days.

# AXIS 500

The Universal Table for All Surgical Discipline

## BASIC POSITIONS



## GYNAECOLOGY / UROLOGY



## NEURO & SPINE SURGERY



CARBON SPINE FRAME "CSF"

SCHURE SPINE FRAME

NEURO HEADREST SYSTEM

## MODULAR TRACTION SYSTEM "MTS"



GENERAL SURGERY

TRACTION FOR FEMUR - 1 LEG

TRACTION FOR FEMUR - 2 LEGS

SHOULDER SURGERY


MIS-BAR - DAA & HIP ARTHROSCOPIES

TRACTION FOR TIBIA

MTS TROLLEY

## Free to Control

The intuitive hand control ergonomically shaped, clearly arranged keypad and LCD display of all important functions.



# SURGI<sup>IRIS</sup>

## X3<sup>MT</sup> LED SURGICAL LIGHTS




Innovative Patented Technology

**FOCUSMATIC**

**X3<sup>MT</sup>** is a LED surgical light consisting of an aluminum open body with three petals.

**LED petal:** the LED sets allow constant central lighting level – regardless of the illuminated field diameter – and colour temperature, while providing a CRI of 98.

**Double yoke or frontal yoke**



**Control handle**

**Sterile central handle:** allows handling and orientation of the cupola by the surgeon.

# SURGIRIS



## EPURE LED SURGICAL LIGHTS



EPURE lights are truly innovative LED surgical lights. Designed and manufactured by Surgiris, they feature the exclusive and patented **FOCUSMATIC** and **MOTION-MATIC** technologies as well as intelligent LED management and control.

Compact and reliable, EPURE surgical lights offer **excellent value for money, making innovation accessible to all.**



Tel Nos.: +63 (2) 8636 3580 | +63 (2) 8637 0930  
Mobiles: +63 917 851 3948 | +63 919 082 6498  
Email: [helpdesk@wellnessproinc.com](mailto:helpdesk@wellnessproinc.com)

#56 San Rafael St., Brgy. Kapitolyo,  
Pasig City 1603 Philippines

[www.wellnessproinc.com](http://www.wellnessproinc.com)



*“Strength comes from the inside.”*

### GL5 Mobile lifter

The all-around safe, reliable & versatile mobile patient lifter with electric width adjustment patient.

Making transfer easier and ergonomically sound, reducing the risk of back injuries to care providers, and improving the comfort and safety for both the carer and user.



### GLS5 Active lifter

An easy, comfortable and safe way to raise and support users with partial or complete control of the upper body



Lifting to and from the floor sitting to lying transfer - vice versa



Lifting to and from the bed Supine transfer/ Horizontal lifting



Heavy lifting / Bariatric care

## GH3 Ceiling Hoist System

The GH3 system is a modular design concept whose classic yet contemporary appearance makes sure GH3 installations blend perfectly with both new and existing buildings, architectural styles and interior design.



### GH3 Lifting Module

Lifting module. User-friendly and extremely powerful and fast lifting module for permanent installation. Lifting capacity 200 or 250kg



### Rail Systems

The rails are what the ceiling hoist is installed in, so they are important in making lifting and moving systems work well. That's why Guldmann rail systems are designed to ensure sure that the hoist moves quietly, smoothly and with a minimum of friction, using as little power as possible. And this smoothness applies right across the capacity range, from 175-500 kg.



### GH1 Ceiling Hoist System

The GH1 ceiling hoist system meets all the day to day lifting requirements encountered in smaller environments such as private homes, group homes or homeless shelters. The discreet, classic design of the GH1 system means it take up very little space.



### GH1 F Lifting Module

With the GH1 F you can decide to install rails in more or all of the rooms and settle with a few lifting modules that in the day-to-day care can easily be moved from one room to another.

# Slings

The comprehensive Guldmann range of slings features slings for many different care tasks, ranging from lifting, positioning and moving to training, mobilisation and activation. The right sling ensures comfort and safety for the user, along with a proper working position for the carer.



## ABC Slings

ABC lifting slings are available in three categories:



**Active**  
Lifting sling that requires active involvement from both carer and user



**Basic**  
Lifting sling that suits most users and most situations



**Custom**  
lifting sling that caters for special needs.

## ABC Kids

ABC lifting slings for children are available in three categories:



**Active**  
children's sling that requires active involvement from both carer and child



**Basic**  
children's sling that suits most children and most situations



**Custom**  
children's sling that caters for the child's special needs

## Other Sling Types Available

Includes Bariatric Slings, Disposable Slings, Turning & Handling



**Bariatric Slings**  
developed with a focus on respect, dignity and comfort, and to meet the challenges that can arise in lifting




**Disposable Slings**  
The disposable slings are delivered in boxes of ten, and are available in both adult and child sizes



**Turning & Handling**  
Handling, positioning and turning slings make a wide range of care procedures easier.

# COLLAPSED LUNG?

## Drain Easy with SINAPI



### 1 Product Safety

- No clamping needed
- No backflow
- Built-in suction regulator for both high positive pressure valve and high negativity vent
- No exposure to bodily fluids at any time
- Provides feedback on actual suction strength

### 2 Drainage Efficiency

- Minimal dead space (60ml) within collection chamber for faster regeneration of negative intra-pleural pressure during the breathing cycle = earlier catheter removal.
- Very low opening pressure (1cmH2O) of valve for low resistance to flow despite the amount of fluids drained = high efficiency of air or fluid drainage from the chest.
- Confirmed suction action via in-line manually depressed suction bulb = safe disconnection of patient from vacuum source for early mobilization.

### 3 Early Patient Mobilization

- Light and compact construction for easy carry
- Orientable in any direction
- Wearable via shoulder sling shot

### 4 Easy Operation

- Easy to set up
- Allows for easy milking of tube
- No fear of backflow
- Drain via built-in tap
- Built-in bed hooks (for some models)
- Little maneuvers needed to monitor and manage

### 5 Cost-saving

- Only one unit for patient's hospital stay
- Efficiently evacuates fluids/air for shorter hospital stay
- Closed system with amenities for less patient complications
- Space saver in medical waste containers

#### FEATURES:



Drainage bag with bayonet connector



#### KEY BENEFITS

1. Product Safety
2. Drainage Efficiency
3. Early Patient Mobilization
4. Easy Operation
5. Cost-Saving

#### KEY SPECS

- Silicone Material Tubing
- Schaeffler Valve
- Ball Movement for Air Leak Intensity
- 25ml or less water to use
- Compact Construction

## PRODUCT RANGE

**ACCESSORIES:** Stand • Shoulder Sling • Y-Connector • Drainage Bag (D-001)



XS50 - Outpatient



XL200S - Pediatric Suction Control



XL200SC - Pediatric Suction Control & HNV



XL1000SC - Cardiac Suction Control & HNV



XL2000S - Cardiac Suction Control & HNV



XL2000SD - Cardiac Suction Control & HNV



D001 - Drainage Bag

# Bidop™ ES-100V3

*Bi-Directional Doppler  
Hand Held Doppler with  
waveform LCD*

- 30 memory storage
- USB computer interface & optional software
- Optional PPG (photoplethysmography)
- Auto gain control



# Mini Doppler ES-100VX

*Pocket Doppler  
with Great Sounds*

- Pocket Doppler with 40% improved GREAT SOUNDS
- Improved volume control for easy operation
- 9V alkaline battery



Elegance  
Meets Flexibility

# NEPTUNE

## Touchless Dispenser



**Reliable & Hassle Free Touchless  
Innovation for Effective Infection Control**

(with Breakthrough OPHARDT Technology at its core)

Ivory

### COLOR VARIANTS



White



Black

\*Standard units come with plain windows



Wide range of sizes



ADA compliant



Secure trick lock  
or key



Hygienic disposable  
pump and cartridge  
system



Includes a wall mounting  
plate with tab feature  
for easy removal



High-end photo  
quality label options


**J** Made in  
**Japan**

**A&D**  
A&D Medical  
Tokyo, Japan

## Multi-User Blood Pressure Monitor TM-2657P



- One button, fully automatic operation
- Displays Systolic, Diastolic and Pulse in a large LED screen
- Quick print out of result in 2.3" width thermal paper
- Power Source: 100-240v, 50-60 Hz AC adapter
- Clinically validated accuracy



Accurate Reading

Easy and Fast Operation

Result Printout

Sanitary Solution

Customization and Integration

perfect fit for installation in clinics, assisted living centers, work place, or many other environments where multiple users need to get fast, comfortable, and accurate blood pressure readings.

# UM-102

## Mercury-Free Sphygmomanometer

- Digital readout for pulse rate
- Cuff holder and grip for carrying
- Durable and chemical resistant body
- Full ABS housing light weight body
- Anti-glare display minimizes reflected light



## UM-201

### Simple Operation BP monitor

- Durable and chemical resistant body
- One-touch operation
- Large LCD display
- Irregular Heartbeat (IHB) Indicator



## UM-211

### All-in-one BP monitor

- Dual measurement mode (Auscultatory & Oscillometric)
- Room temperature display
- Cuff holder & grip for carrying

**Made in Japan**

**AND**  
A&D Medical  
Tokyo, Japan

# UA-767JP

## Digital Blood Pressure Monitor

- Clinically validated accuracy
- Robust, long time use
- Whole family use
- Recommended by Filipino MD's



Clinically Validated



Original AC Adapter Included



IHB Indicator



Cuff Fit Error Indicator



Memory

WHO

WHO Blood Pressure Classification Indicator



SlimFit Cuff/Wide Range Cuff



Movement Error Indicator



Clinically Validated

# UA-704

## PalmTop Blood Pressure Monitor

- Innovative "PalmTop" Style
- Irregular Heart Beat (IHB) indicator
- SlimFit Comfortable Cuff
- Uses Single AA Battery



IHB Indicator



SlimFit Cuff/Wide Range Cuff



Weighs only 73 grams



Uses single AA Battery

## Trustworthy Blood Pressure Readings



Tel Nos.: +63 (2) 8636 3580 | +63 (2) 8637 0930  
Mobiles: +63 917 851 3948 | +63 919 082 6498  
Email: helpdesk@wellnessproinc.com

#56 San Rafael St., Brgy. Kapitolyo,  
Pasig City 1603 Philippines

www.wellnessproinc.com

# Baby Shuttle

## Transport Incubator

3 Variants Available:

Normal

Special

Intensive Care

Alarms for

High & Low Temperature

Low Battery

Max Temperature

Fan Failure

Main Failure

- Light in weight (fully loaded at 50kg)
- Self-loading trolley (Good for Helicopter, Plane, Ambulance)
- 2 Ways to access baby - by his head or by his side
- Stand alone Battery of 3.5 hours



Servo Control

Skin Temperature

Air Temperature



# Polytrend

## Intensive Care Incubator

Next Level Infection Control Features

Round Hood Corners

Doors can open upwards

Monitors

Thermoregulation of the newborn

Air Temperature

Skin Temperature


Peripheral Temperature

Humidity Level

Oxygen concentration level

SpO2 heart rate (optional)

- Largest Internal Hood Space Available
- Completely Autoclavable for Secured Infection Control (Hood, Bed, Body)
- Upgradeable on Plug & Play Basis
- Break Resistant Polycarbonate for Hood/Bed
- Gentle Air Buddle Around Baby
- Ginevri Single Wall Structure & 118 Built in Air Vents
- Comes with Trolley (w/ or w/out height adjustment)
- Variants Available with different options



# Alhena Plus

## Infant Warmer

3 Waves to Control Heating Level

Automatic Mode

Manual mode

Pre-Heating Mode

- 360° Rotatable / Tilttable bed for optimized maneuvering
- Power LEDs are highly directional, irradiating only the treatment area
- Collapsible sides
- Lifespan of lamps long 20,000 hours
- Versatile w/ Lamp Types Built-in
- Skin Temperature Sensor to ensure "Safe for Baby" heat control



# Naos Plus

## Phototherapy Lamp

- Equipped with 6 Power Blue LEDs
- Power LEDs for treatment of hyperbilirubins in newborns
- Manages multiple therapy cycles



# Mira

## Fiber Optic Phototherapy Lamp

- Equipped with 6 Power Blue LEDs
- Power LEDs for treatment of hyperbilirubins in newborns
- Manages multiple therapy cycles

# Baby Start System

Infant T-Piece Resuscitator



*“Crucial on hand in first 3 minutes of premature baby’s arrival”*

## Features

- Simple and efficient device to resuscitate asphyxiated neonates by mask or endotracheal tube.
- Resuscitation pressures (PIP and PEEP) are easily readable on the fast response manometer.
- Uses a single use patient circuit completed with an expiratory valve and silicone neonatal mask.
- Can be installed, by dedicated support, on Infant Warmers, on Transport Incubator or Resuscitation Trolley.

## Key Components of System

Baby Start



WETTY Humidifier



Oxygen Gas Blender



Single Flowmeter



# MASTERMED® C

PREMIUM BLOOD PRESSURE MEASURING DEVICES



The KaWe MASTERMED® C Premium Blood Pressure Monitoring provide clinicians the highest level of quality and value.

The series is defined by its innovative and practical features, smooth surfaces for easy cleaning, as well as its particularly precise gauge. These blood pressure monitoring devices are available as wall-mounted, table-mounted and freestanding models.

These units have a double-walled housing that ensures additional stability and protection and are therefore designed to withstand twice the pressure. Each model comes with a generously-sized cuff holder.

## AVAILABLE MODEL



### WALL-MOUNTED MODEL

#### Features:

- Swivelling mounting bracket
- For use with the wall-mounted KaWe MedCenter 5000 work station
- Screws and wall anchors included
- Spiral hose with a range of 3 meters



### TABLE TOP MODEL

#### Features:

- Integrated handle and bases that prevents scratches
- Spiral hose with a range of 1.2 meters
- Cuff, 2 tube



### STAND MODEL

#### Features:

- Stable 5-foot roller stand made of high-quality stainless and sturdy plastic juncture
- Adjustable height: 85 to 130 cm
- Quiet Ø50 mm double castors: 2 locking castors, 2 safety castors and one anti-static castor
- Spiral hose with a range of 3 meters

## FEATURES

- Complies with ISO 81060-1 standards
- Clear gauge with thin, red dial for exact measurement readings
- Numbered instrument dial reduces the risk of parallax
- 150 mm dial diameter
- Scaling: 0 - 300 mmHg
- Beryllium-copper diaphragm maintains the calibration up to 600 mmHg
- Swivel-mounted
- Equipped with PU/nylon-Velcro cuff for adults
- Plug-in adapter for quicker cuff exchange
- Stays calibrated for 5 years, guaranteed
- 2-year guarantee on device when used properly

# OTOSCOPE

Piccolight Fiber Optic LED



night.  
sky.  
stone.

## FEATURES

- Fiber Optic LED illumination
- Illuminance approx. 8,000 Lux
- Lifetime bulb approx. 100,000 h
- Shockproof and energy-saving
- Battery handle AA
- Disposable ear funnels (24)  
10 x Ø 2.5 mm, 10 x Ø 4.0 mm
- In cloth bag
- Battery: 2 Mignon (type AA)

# DIAGNOSTIC SET

## BASIC SET C10 / E10

OTOSCOPE SET /  
OPHTHALMOSCOPE SET



### SET INCLUSIONS

Otoscope Head

Ophthalmoscope Head

Handle

### ACCESSORIES

- 3 usable ear funnels
- 2 laryngeal mirrors
- Curved lamp support
- Blade holder made of plastic
- Nasal speculum, chrome plated  
with set screw
- Spare bulb
- Packaged in a white hard-top box

MEDIZINTECHNIK  
seit 1890

**KaWe**  
GERMANY

# TOP CARDIOLOGY

## FLAT-HEAD STETHOSCOPE



Name Plate



Ear Tips



Available in colors of

■ **Burgundy** ■ **Black**

### FEATURES

- Flat - head stethoscope
- Stainless steel
- Diaphragm Ø45 mm
- Total length approx. 64 cm
- Single-tube system with integrated head piece made of stainless steel
- The ideal stethoscope for cardiologists
- Includes spare parts set
- Weight: approx. 210 grams



Item No.  
701145



# ERO-SCAN® Monochrome

## OAE Test System


- Performing screening and diagnostic testing with TE and / or DPOAE
- Fast automatic test with Pass / Refer and graphical test result
- Sharp, colored organic LED Display
- Direct evaluation via value and bar diagram
- Up to 12 frequencies displayed within the device, DPOAE up to 12kHz
- High noise immunity for operation in normal clinical environment
- Lightweight, small earprobe
- Environmentally friendly due to long life, rechargeable battery
- Bluetooth® communication to optional printer and PC
- Evaluation of results in PC
- Multiple, easily selected language options

# MAICO MB 11 BERAphone®

## AABR Test System

- Automatic, Reliable ABR-screening Results
- Innovative BERAphone design given cost saving integrated electrodes
- CE-Chirp-Stimulus / Impedance Check for fast results
- Automatic Impedance Check indicates impedance conditions
- PC Connectivity for Power / Data Transfer

User-friendly  
software



Item No.  
701174

The special evaluation software performs the test automatically and the Pass or Refer Result will show within a few seconds.



## VITROMAP

*Digital Desktop Pulse Oximeter*

The VITROmap is a budget-priced pulse oximeter that works up to 12 hours on the integrated, rechargeable battery or directly on a mains connection.

- Lightweight, Compact and Attractive Design
- Adjustable Alarm Limits
- Simple Operation



## PEARL 8

*Digital Handheld Pulse Oximeter*

The PEARL8 Pulse Oximeter with its new digital signal processing technology works reliably and can be used from neonate to adult.

- Reliable measurement of SpO2 and pulse rate
- High quality and brilliant contrast TFT screen
- Can be used on adults, children and neonates
- Patient Monitoring



## NANOXECO

*Digital Handheld Pulse Oximeter*

The NANOXeco pulse oximeter has everything you expect to see in a handheld device used for spot checks, at an extremely attractive price.

- Measures a patient's SpO2 and pulse rate
- Large, bright and easily readable LED displays
- SpO2 dependent pulse tone (muteable)

### Accessories



Fingerclip  
Probe P-200



Finger  
Probe R-200



Finger  
Probe RS-200



Small Finger  
Probe PS-200



Universal  
Probe Y-200



Wrap  
Probe WR-200

# Capnograph

CO<sub>2</sub> Monitored Accurately

medlab  
medizinische Diagnosegeräte GmbH



## CAP10

*Sidestream Capnograph*

The new CAP10 Capnograph with sidestream technology fulfills all requirements you have in clinical and pre-clinical applications for monitoring a patient's breath and CO<sub>2</sub> status.

- Lightweight, compact and attractive design
- Sidestream system with integrated pneumatics
- High resolution and bright colour TFT screen



## CAPNOS/CAPNOX

*Digital Handheld Capnograph*

The CAPNOS is a compact, mobile capnograph for ambulances, day-surgeries and clinical environments.

The unit continuously monitors the expiratory carbon dioxide concentration and the respiration rate of a patient

Measurement of the etCO<sub>2</sub> values in a patient's breath

- Display of etCO<sub>2</sub>, Respiration Rate and Capnogram
- Trend Display of two Parameters
- Operates w/ a rechargeable battery pack 10-12 hours
- Optional with SPO<sub>2</sub> function (CAPNOX )

### Capnography

is the monitoring of the concentration or partial pressure of carbon dioxide (CO<sub>2</sub>) in the respiratory gases. Its main development has been as a monitoring tool for use during anesthesia and intensive care. It is usually presented as a graph of expiratory CO<sub>2</sub>.



# seca

Precision for health

## Medical Body Composition Analyzer

# mBCA

medical Body Composition Analyzer

Simplified work, quick measurement  
and medically relevant results



## “Clinic-level results”

### The seca mBCA

- Standing aid
- Automatic electrode detection
- Flat platform
- Short measurement time
- Easy entry of patient data
- Flexible storage concept
- Simple user navigation
- Integration in every EMR system
- seca 360° wireless

## seca mBCA 515



# seca mBCA 525

Results in 30 seconds!

Medical Body Composition Analyzer  
for determining body composition whilst in a supine position

*Measurement of fat, water, muscle and other body parameters in a supine position*

*Perfect for mobile use at only 3 kg with long-lasting rechargeable battery*

*WiFi connection between measuring mat and monitor guarantees reliable data transmission*

*Easy-to-understand graphic presentation of measurement results on the touchscreen monitor*

*Medically precise, validated against reference methods of the gold standard*

*Internal storage for results of approximately 100,000 measurements*



“For supine or seated patients”



### Lightweight on the go.

Its mobility is optimized even further through the mobile stand.



### Perfectly displayed results.

The clearly organized presentation of all parameters makes diagnosis extremely easy and simplifies the patient consultation.



### Always ready for use.

The measuring mat is charged by induction in the practical transport compartment.




### Reliable measurement.

The correct placement of the electrodes is illustrated on the measuring mat.

### Easy and quick use.

The Easy Snap™-electrodes can be applied very quickly and the square adhesive electrodes also simplify the correct placement.



## Multifunctional Wheelchair Scales

# seca

Precision for health

### seca 684

Wireless multifunctional scale with stable handrail and high load bearing capacity with castors for ease of transport

*Patient can stand or sit when being weighed*

*The rotating display can be oriented towards the patient or medical personnel*



### seca 656

Wireless electronic platform scale with innovative memory function

*Especially spacious platform for all common roll-in systems*

*Stable construction with very high capacity*

*Transportable due to castors and integrated handles*

*Innovative memory function for determining precise weight of patient*



*Precision weighing in all situations*



## Multifunctional Wheelchair Scales

**seca**  
Precision for health

### seca 676

Electronic Wheelchair Scales  
with Handrail & Transport Castors

*Integrated ramp*

*Castors & Folding feature for easy transport*

*360 kg weighing capacity*

*Precise weighing even with  
handrail support (innovatively SECA)*



### seca 664

Digital wheelchair / Platform Scales

*Swiveling LCD with large letters for easy reading*

*No side rails allowing for placement  
in more (possibly tight) places*

*360 kg capacity*

*Value for Money Digital wheelchair / platform scales*



# Wireless Ultrasonic Measuring Station for height, weight and BMI

**seca**  
Precision for health

*“ With voice guidance.”*

## seca 286

Wireless ultrasonic measuring station for height and weight with voice guidance

*Fully automated measurement is possible through the optional voice guidance*

*Three pairs of sensors and constant auto-calibration ensure the most precise results*

*The unrivaled digital evaluation algorithm eliminates imprecision for the highest measurement accuracy*

*A wide variety of output options for the measurement results*



### Sensors

Three unique pairs of sensors transmit and receive ultrasonic signals separately and thereby ensuring the most precise measurement results.



### Display

The measured values on the display can be switched on or off and has a keylock.




### Wireless printer

A wireless printer enables the results to be easily taken to the nurse or physician directly after measuring.

### Non-slip tread

The scale has a non-slip tread made of easy-to-clean safety glass.



# Measuring Stations and Column Scales

# seca

Precision for health



## seca 786

Mechanical Column Scales with Large Round Dial

*Effortless weighing*

*Determine Patient Weight As Soon As He Steps On the Platform*

*with height rod*



## seca 700X

Mechanical Column Scale with Eye Level Beam & Attached **Seca 224** Height Rod

*220 kg Capacity (50 g increment) / 60-200 cm range*

*Non Fading / Engraved Numbers Assuring Always Easy-to-read of Results*

*SECA Breakthrough Innovations on the Familiar Beam Scale*

*The broad and flat construction provides easy step-on access and a firm foothold*



## seca 755

Mechanical Column Scales with large dial and automatic display of the BMI range in kg or lbs version

*Quick and easy ascertainment of BMI*

*Simultaneous evaluation of nutritional condition*

*Clearly designed dial and low platform*



# Body Scales

# seca

Precision for health

## seca 374x

Wireless baby scale w/ shell-shaped tray, raised LCD and Measuring Rod



*Wireless Communication  
Best for Hospital Nurseries & Pediatrician Clinics*

## seca 727

Wireless baby scale with fine 1g graduation



*Best for NICU*

*For weighing babies in a sling seat suspended under the scales*



## seca 310

Mechanical circular dial scale



*Removable Tray*



## seca 354 Cap. 20kg

Multi-purpose, electronic baby and infant scale

## Body Scales

### seca 417

Light and stable measuring board for mobile use



Foldable, Stable, Light Mobile  
Measuring Board for Babies  
with 1mm increment

*Used & Deployed by UNICEF  
Worldwide*




### seca 210

Mobile, Rollable, Washable,  
Non-Fading Measuring Mat  
for Babies with 5mm increment

### seca 416

Infantometer for Babies (up to  
24mo.) with locking system



### seca 212

Measuring Tape for Head  
Circumference of Babies  
and Toddlers

*Precise and easy-to-read  
measurements*



### seca 207

Measuring rod with extra wide measuring slides

*measuring rod can be firmly screwed to the changing table  
0 - 99 cm measuring range*

## Flat Scales

**seca**  
Precision for health

### seca 874

Electronic flat scale for mobile use with push buttons and double display

*Extremely robust construction*

*Easy-to-use design with long-life batteries*

*Precision unaffected by humidity levels*

*Double display lets the patient and medical personnel read the results from two different prospects at the same time*



Both doctor and patient can see the result at the same time in the double display



### seca 803

Electronic flat scale with very high bearing capacity of up to 200 kg



### seca 760

Robust Mechanical Dial Floor Scale  
*Also available in White, Blue & Yellow*



### seca 813

Electronic flat scale with very high bearing capacity of up to 200 kg

### seca 856

Organ & Diaper Scale  
*Capacity: 5 kg*



### seca 852

Digital portion & diet scale  
*Capacity: 3 kg*





## seca 206

Microtoise / Mechanical Measuring Tape with Anti-Twisting Construction for longer lasting use

*Easy to read height  
Durable metal tape*

## seca 201

Ergonomic circumference measuring tape



## seca 213

Portable Stadiometer



*only 2.4 kg, the seca 213 mobile stadiometer is especially suitable for use in the field but also ideal for doctors' practices and pediatric hospitals*



*Products designed down to the last detail*

# Electric Safety Testers

## SECULIFE | ST w/ (P)Si Module M6930 Test Instrument

As a multifunctional test instrument, the Seculife ST is ideally suited for performing and documenting measurements at medical devices. It allows for precision testing of portable medical devices, and can also be used for large devices such as X-ray equipment and MRI scanners.

The SECULIFE | ST detects mains supply power, safety class and any mains connection errors at the device under test. The test sequence is significantly simplified by means of the automatic mode :

- Select a regulation.
- Press the start key.
- Read the results.

### Features

- Mains Plug Polarity Reversal
- Automatic Recognition
- Display
- Menu Driven Test Sequences
- Expandable
- Basic Instrument and Expansion Features
- Leakage Current Test



Easy-to-operate

The SECULIFE ST is the only measuring instrument available on the market which has been designed to test medical devices operated with 3-phase electrical current. The three-phase current adapters AT3-IIIE and AT3-IIS 3-phase current adapters can be furnished to this end as optional accessories.

#### DAkkS Calibration Certificate included

Power shutdown as soon as the leakage current exceeds approx. 15mA provides a maximum of safety for the user

High-voltage test for:  
DIN EN 60950, DIN EN 61010, DIN EN 60335 and DIN EN 60601



German  
Accreditation Body  
D-K-15080-01-01

DAkkS Calibration Certificate as Standard Feature

## SECULIFE | HIT MD

TRMS System Multimeter for Use in Medical and Hygienic Sensible Ranges

*The measurement and testing of measuring equipment in the field of medical technology and food and pharma industries necessitates absolute reliability and a broad range of applications.*



IP65



German Accreditation Body  
D-K-15080-01-01

DAKKS Calibration Certificate as Standard Feature

## Features

- 23 multimeter functions
- Voltage measurement
- Auto-ranging current measurement from 100µA
- Capacitance and resistance measurement
- Measuring categories: 600V CAT III and 300V CAT IV
- 1 kHz low-pass filter
- TRMS AC and AC+DC, 20 kHz bandwidth
- Data storage for more than 15,000 measured values
- Hygiene certificate from University of Marburg

### • Interfaces:

**Bidirectional infrared interface for communication**  
*Optionally available IR-USB adapter*

### • Power Supply:

**Battery operation**  
**Mains operation via optional broad range variable power pack**



**SECULIFE | SR**  
Test Instrument

modular test system is the ideal solution for complying with medical device legislation: diverse, individual and cost-saving.



**SECULIFE | SB**  
Safety Analyzer for Medical Devices with Integrated Database

offers all of the advantages of a universal, computer controllable measuring instrument for the field of medical technology



**SECULIFE | IP**  
Measuring Instrument

instrument is capable of executing all measurements for testing the effectiveness of safety measures in electrical systems

# High-Frequency Surgery Measuring Instruments

## SECULIFE | ESPRIME

Multimode Electrosurgical Analyzer

The new SECULIFE ESPRIME is in a class of its own, and its real time operating system offers virtually unlimited expansion capabilities in the future that make it the ONLY logical choice for customers looking for a full-featured analyzer.

### Features

- Capture, Store, Print RF Waveforms
- Industry Standard RF Current Measurement for multiple pulsed waveforms
- Automated User-Definable Testing Sequences
- Ultra-fast digitization of the complex RF waveform
- USB- (3), RS232-, and Ethernet communication Ports
- Remote communications capability with ESU Generators
- Embedded Real-Time operating System with 1/4 VGA color touch screen display



## SECULIFE | ESXTRA

Electrosurgical Analyzer



## SECULIFE | ES PRO

Electrosurgical Analyzer

is an intelligent instrument for testing the latest generation of defibrillators



## SECULIFE | ES BASE

Multimode Defibrillator Analyzer

is an easy-to-use power meter for testing high frequency generators used in the field of HF surgery

# Photometry Instruments

## SECULIFE | IT-B | IT-C | IT-BASE

Luxmeter

The SECULIFE IT BASE is a high precision luxmeter for use in hospitals and medical practices. It reliably measures the illuminance of daylight, as well as all types of artificial light sources, and is ideal for quality assurance: precision in accordance with class C makes it suitable for a broad range of measuring tasks.



### Features

- Accuracy amounts to  $\pm 3\% \pm 1$  digit display value
- Data memory for up to 100 single measuring values
- Illuminance can be measured in Lux or Footcandle
- Wide measuring ranges
- Easy operation
- High resolution

### SECULIFE | IA

Luminance Analyzer

light measurement technology for medical reporting



### SECULIFE | IM

Luminance Analyzer

for measuring luminance



### SECULIFE | IS

Luminance Analyzer

Illuminance Sensor for the Ambient Light Control in Reading Rooms



# Function Testers

## SECULIFE | DF BASE

Defibrillator Analyzer



The SECULIFE | DFBASE is an easy-to-use power meter for testing high frequency generators used in the field of HF surgery. Measured values are indicated at a large analog display. Works without any external power supply and without batteries.

## SECULIFE | DF PRO

Multimode Defibrillator Analyzer



The SECULIFE | DFPRO multimode defibrillator analyzer is an intelligent instrument for testing the latest generation of defibrillators. It measures discharged energy and provides information regarding the defibrillator pulse.

# Patient Simulators

## SECULIFE | BP PRO

NIBP Simulator for Testing Patient Monitors

The SECULIFE | BP PRO is the high-end product from the new range of NIBP simulators and is distinguished by a great variety of functions. It includes functions such as NIBP, IBP, ECG, temperature, arrhythmia, respiration, leak detection and much more.

### Features

- Optional peak pressure detect with simple reset
- Synchronized invasive blood pressure output
- ±1% of reading pressure accuracy
- ECG output with full NSR waveform
- Digital pressure envelope offset
- With multifunctional connection
- Small, handy, lightweight
- Flash programmable



### SECULIFE | NIBP

Multi-parameter Simulator

for testing patient monitors



### SECULIFE | PS300

Multi-patient Simulator for ECG

for testing patient monitors



### SECULIFE | PS200

Multi-patient Simulator for ECG

for testing patient monitors



### SECULIFE | PS100

Multi-patient Simulator for ECG

for testing patient monitors



# Measuring Instruments



## SECULIFE | IF+ M695D

Infusion Pump Analyzer

The Seculife IF+ is designed to test the flow rates of intravenous infusion pumps to ensure their correct operation. It is able to test up to two pumps simultaneously.

All tests are controlled by a Microprocessor which calculates and displays the results. The fluid path is unobstructed, making them easy to clean. The tubes are plastic, not glass so they are rugged. The cables and chambers are interchangeable.

## Features

- Simultaneous operation of two chambers
- $\pm 1\%$  of reading flow accuracy
- Digital calibration – no pots to turn
- Large graphics display with cursor selection
- Easy to use and clean
- Auto chamber recognition



**SECULIFE | DP PRO**  
Precision Digital Measuring Instrument

is a high-precision digital pressure gauge and measures compatible gas and liquid pressure values in various units of measures



**SECULIFE | DP BASE**  
Precision Digital Measuring Instrument

is a microprocessor controlled digital pressure gauge which permits positive as well as negative pressure measurement in liquids and gases using various physical units of measure



## Type FM8000G

### FEATURES

- Light weight, easy operation
- 8.4" screen color LCD display, rotatable screen to 60
- FHR 50-240 bpm
- Manual/Automatic records fetal movement
- Sound and color alarm for high and low fetal heart rate
- Continuous 24-hour real-time monitoring optional
- Continuous 12-hour patient curve and date storage, playback print
- With picture freeze function
- Single, Twins Monitoring optional
- Extra-long life with built thermal recorder



## Type FM8010

### FEATURES

- Color TFT folding color LCD : 0.2" color TFT folding LCD
- 9 crystal pulse wave Doppler transducer
- 152mm wide paper (US/International standard)
- Built-in rechargeable battery
- Accurate fetal movement detection and accounting
- Central monitoring interface

## OTHER PRODUCTS

### ECG MONITOR



**CARDIPIA 400**  
Type 406N/406R

- Compact and light weight
- 4-lead patient cable
- Clear and reliable ECG trace
- Built-in rechargeable battery



**CARDIPIA 200**  
Type 201N/203N

- Compact and light weight
- 4-lead patient cable
- Clear and reliable ECG trace
- Built-in rechargeable battery



**A&D**  
A&D Medical  
Digital Blood Pressure Monitor

**Made in Japan**

Use Quality Assured  
Medical Device Only.

 Clinically validated as A/A  
by European Society of Hypertension

 with **FREE**  
Original AC Adapter




**NAVEH PHARMA**

**CleanEars**  
For Ear Wax Removal  
Wash & Dissolve Jet Spray

**CleanEars**  
For Ear Wax Removal  
Wash & Dissolve Jet Spray

- PREVENTS EAR PLUGGING
- NO STEROIDS
- SAFE AND EFFECTIVE

**CleanEars**  
For Ear Wax Removal  
Wash & Dissolve Jet Spray

Get yours at

 **mercury drug**  
NAKASISIGURO GAMOT AY LAGING BAGO

 **Lazada**

Exclusive Philippine Distributor  
 **WELLNESSPRO**  
☎ (02) 8636-3580 | www.wellnessproinc.com